

Edio Mauricio Llerena Mancías - A00828309  
Jesus Gabriel Leos Perales - A01412140  
Andres Martin Vivanco Palacios - A01705733  
Francisco Valadez Rojas - A01610820

### MRI activity.

Open the simulator in <https://phet.colorado.edu/en/simulations/mri>

#### 1. NMR

Place yourself in the tab *Simplified NMR*. Remember that the Larmour relationship relates the magnetic field to the resonant frequency:

$$\nu_L = \frac{\gamma}{2\pi} B_0$$

where  $\nu_L$  is the resonant frequency,  $\frac{\gamma}{2\pi}$  is the gyromagnetic ratio and  $B_0$  is the magnetic field.

Table1. Gyromagnetic ratios for different nuclei.

Nuclei	Gyromagnetic ratio	Nuclei	Gyromagnetic ratio
$^1\text{H}$	42,58	$^{65}\text{Cu}$	12,09
$^7\text{Li}$	16,55	$^{75}\text{As}$	7,291
$^9\text{Be}$	5,984	$^{77}\text{Se}$	8,118
$^{11}\text{B}$	13,66	$^{81}\text{Br}$	11,50
$^{13}\text{C}$	10,71	$^{87}\text{Rb}$	13,93
$^{15}\text{N}$	4,314	$^{93}\text{Nb}$	10,41

<sup>17</sup> O	5,772	<sup>117</sup> Sn	15,17
<sup>19</sup> F	40.05	<sup>121</sup> Sb	10,19
<sup>23</sup> Na	11,42	<sup>127</sup> I	8,518
<sup>27</sup> Al	11,09	<sup>133</sup> Cs	5,584
<sup>29</sup> Si	8,458	<sup>195</sup> Pt	9,153
<sup>31</sup> P	17,24	<sup>199</sup> Hg	7,590
<sup>35</sup> Cl	4,172	<sup>203</sup> Tl	24,33
<sup>51</sup> V	11,19	<sup>207</sup> Pb	8,907
<sup>55</sup> Mn	10,50	<sup>209</sup> Bi	6,841
<sup>59</sup> Co	10,05		

Use the Larmour relationship and the gyromagnetic ratios of various nuclei shown in table 1 to complete table 2. Check your results in the simulation by setting the appropriate frequencies and main magnet field, take a screenshot of the nuclei emitting energy to include in the report. Try to find the last nuclei (????) by playing with the simulation and register the frequency at two different magnetic fields.

Table 2. Different settings to achieve energy emission.

Nuclei	Magnetic Field	Resonant Frequency	Magnetic Field	Resonant Frequency
Hydrogen	0.75	31.935	2.5	106.45
Nitrogen	2.5	10.785	-	-
Sodium	1.515	17.3	2.75	31.405
Carbon-13	1.75	18.7425	2.5	26.775
Oxygen	1.999	11.54	3.0	14.316
Sulfur	3.0	~10	-	-
Copper ( <sup>65</sup> Cu)	3.0	36.3		

Simplified MRI (1.08)

File Help

Simplified NMR Simplified MRI

Legend

- NMR Nuclear magnetic resonance
- MRI Magnetic resonance imaging
- Atomic nucleus
- Magnetic Field
- Radio Wave

Energy

Main Magnet

Field: 0.75 Tesla

Sample Atom

hydrogen

Power

0 25 50 75 100

51 %

Radiowave Source

Wave view Photon view


Frequency


10 35 60 85 110


31.935 MHz



Pause Play

The simulation interface displays a central grid of 20 atomic nuclei, each represented by a red and blue sphere with a small arrow indicating its magnetic moment. The nuclei are arranged in a 4x5 grid. Above and below the grid are yellow bars representing the main magnet. The interface includes a 'Radiowave Source' panel at the bottom with a 'Power' slider set to 51% and a 'Frequency' slider set to 31.935 MHz. A 'Main Magnet' panel on the right shows a 'Field' slider set to 0.75 Tesla. A 'Sample Atom' dropdown menu is also present, currently set to 'hydrogen'. A legend on the right explains the symbols used. An energy level diagram shows two discrete energy levels with nuclei in different states. At the bottom, there are pause and play buttons.



























## 2. MRI

Move to the *Simplified MRI* tab

1. Set the *main magnet field* to 1.0 Tesla, leave the *gradient magnets* in zeros, activate only *show head*, and *show magnetic field* (be sure that *show atomic nuclei* is deactivated), set the frequency in 43 MHz. Finally set the *power* to 50% and observe the flow and distribution of the emissions. After a while observing the emissions, click on *add tumor*, wait for around 7 seconds so the flow distribution stabilizes, look at how the emission changed and try to guess where the tumor is located.



Explain how the emission allowed you to find the correct location: We can identify the tumor by looking at the dark lines (signals) going left to right, the tumor is located where there is a higher concentration of these signals. In the second image we can say that the tumor is located in the middle left side of the patient's head.

2. Play with the main magnet field, frequency, and gradient magnets (both, horizontal and vertical) to try to obtain an emission focused mainly in the zone of the tumor (register your best guess, it doesn't need to be perfect). Answer the following questions.

Best guess: main magnet: 1 Tesla horizontal gradient: 0.05 Tesla vertical gradient: 0.075  
frequency: 38MHz

What happens when the horizontal gradient increases its magnitude? How does it affect the emissions? The magnetic field is altered, causing it to increase as it approaches the detector. This change in the magnetic field, that was originally uniform, needs a change in

the frequency to spin the atomic nuclei and facilitate the emission focus on the tumor we are looking for.

What about vertical gradients? Just as the horizontal gradient provoke alterations on the magnetic field, also the vertical gradient eliminates the uniformity adding external sources that increase the magnetic field from the top to the bottom. These changes allow us to focalize zones and acquire better emissions. If we apply these alterations it is necessary to find the adequate frequency that matches the zone we are looking for.